



UNIVERSITÀ
DEGLI STUDI DI TRIESTE

Dipartimento di Scienze Chimiche e Farmaceutiche

Trieste, 14.11.2022

Descrizione	Quantità
PPL Lined Hydrothermal Autoclave reactors of 25ml	2
PPL Lined Hydrothermal Autoclave reactors of 100 ml	2
PPL liners of 25 ml	1
PPL liners of 100 ml	1

Università degli Studi di Trieste

Dipartimento di Scienze Chimiche e Farmaceutiche
sedi: - via Licio Giorgieri, 1 (amministrazione)
- piazzale Europa, 1
34127 Trieste

dipdscf@units.it
tel +39 040 558 3902
tel +39 040 574181-5587943

www.dscf.units.it
fax +39 040 558 3903
fax +39 040 52572

14.11.2022, 08:58

Почта — KOST VERONIKA [PHD0600114] — Outlook

Quotation_PPL Lined Hydrothermal Autoclave

KOST VERONIKA [PHD0600114] <VERONIKA.KOST@phd.units.it>

08.11.2022, Вт, 11:53

Кому: info@shilpent.com <info@shilpent.com>

Dear Sir/Madam,

Good morning!

My name is Veronika Kost and I am a researcher at the University of Trieste (Italy). I am interested in the following items: PPL Lined Hydrothermal Autoclave (1 reactor+linker of 25 ml, and 1 reactor+linker of 100 ml, 2 reactors in total).

Could I ask you to send the quotation for these items?

With best regards,

Veronika Kost

QUOTATION

Quotation # 8207



Shilpa Enterprises

No. 247, Yadav Nagar, Post Uppalwadi,
Near Automotive square,
440026, Maharashtra, India
+91-9096433131, +91-9168009666
info.shilpent@gmail.com
www.shilpent.com
GSTIN 27AKXPC9094M1Z4

Quotation Date : 08/11/2022

Bill To
University of Trieste
Italy.

VERONIKA.KOST@phd.units.it

#	Item & Description	Quantity	Rate	Amount
1	<p>Hydrothermal Autoclave 25 ml Hydrothermal Autoclave Capacity – 25 ml Product Code – AUTO25X Brand – Shilpent</p> <p>Technical specifications:- Safe Temperature: - 180-200°C Working Pressure: ≤ 3MPa Material: high-quality stainless steel Heating and Cooling Rate: ≤ 5°C/min Inner PTFE (Teflon Liner) Chamber</p> <p>HSN: 84192010</p>	1.00	94.00	94.00
2	<p>Hydrothermal Autoclave 100 ml Capacity – 100 ml Product Code – AUTO100X Brand – Shilpent</p> <p>Technical specifications:- Safe Temperature: - 180-200°C Working Pressure: ≤ 3MPa Material: high-quality stainless steel Heating and Cooling Rate: ≤ 5°C/min Inner PTFE (Teflon Liner) Chamber</p> <p>HSN: 84192010</p>	1.00	189.00	189.00
3	<p>Shipping Charges</p> <p>HSN: 84192010</p>	1.00	39.00	39.00

#	Item & Description	Quantity	Rate	Amount
4	Bank Fees	1.00	70.00	70.00
	HSN: 84192010			
			Sub Total	392.00
			IGST0 (0%)	0.00
			Total	\$392.00

Total In Words: ***Dollars Three Hundred
Ninety-Two***

Bank Details:-

Account Name: Shilpa Enterprises

Account no. 389103000000569

Swift Code - IBKLINBBNGP

IFSC Code: - IBKL0000389

Bank name: - IDBI bank

Branch: - Corporate branch, Nagpur

Terms & Conditions

CLIENT CAN PURCHASE SINGLE ITEM ALSO

Payment Terms: - 100% Advance payment

P. O: Purchase order is required

Payment mode – Wire Transfer (T/T), Credit/Debit Card

Delivery period: Within 15-20 Days after confirmation of order

Delivery Mode: Via Express Air

(M/S. Shilpa Enterprises)



14.11.2022, 08:59

Почта — KOST VERONIKA [PHD0600114] — Outlook

Quotation

KOST VERONIKA [PHD0600114] <VERONIKA.KOST@phd.units.it>

08.11.2022, Вт, 9:47

Кому: Techinstro Team <info@techinstro.com>

Dear Sir/Madam,

Good evening!

My name is Veronika Kost and I am a researcher at the University of Trieste (Italy). I am interested in the following items: *Autoclave setup with liner: PPL Lined Hydrothermal Autoclave (Safe Temp – 260°C), 25ml and 100ml in volume.*

Could I ask you to send me the actual quotation for these items?

With best regards,

Veronika Kost



Techinstro Innovations
 Plot no 463, First Floor, Yadav Nagar
 Nagpur Maharashtra 440026, India
 +91-9765849656, +91-8007799090
 info@techinstro.com
 www.techinstro.com
 GSTIN : 27CDEPM2102E1Z4

QUOTATION

Quote To,
University of Trieste
 Italy
 VERONIKA.KOST@phd.units.it

Quotation No. 1202
 Quotation Date 08/11/2022
 Sales person Website

Sr. No.	Product & Description	HSN/SAC	Quantity	Rate	Amount
1	PPL Lined Hydrothermal Autoclave 25ml Capacity – 25 ml Product Code - Autoclave-PPL-TI025 Brand – TECHINSTRO Technical specifications:- Working Temperature: 280°C Safe Temperature: - 260°C Working Pressure: ≤ 3MPa Heating and Cooling Rate: ≤ 5°C/min Material: Shell made of high quality stainless steel Inner Shell: PPL Liner	84192010	1.00	229.00	229.00
2	PPL Hydrothermal Autoclave 100ml Capacity – 100 ml Product Code - Autoclave-PPL-TI0100 Brand – TECHINSTRO Technical specifications:- Working Temperature: 280°C Safe Temperature: - 260°C Working Pressure: ≤ 3MPa Heating and Cooling Rate: ≤ 5°C/min Material: Shell made of high quality stainless steel Inner Shell: PPL Liner	84192010	1.00	379.00	379.00
3	Shipping Charges	84192010	1.00	77.00	77.00
4	Bank Fees	84192010	1.00	39.00	39.00
				Sub Total	724.00
				IGST0 (0%)	0.00
				Total	\$724.00


Total In Words: **United States Dollar Seven Hundred Twenty-Four**

Bank Account Details:-
 Account Name: **Techinstro Innovations**
 Account no. 40581536347

IFSC Code: - SBIN0000432
Swift Code - SBININBB924
Bank Name: - State Bank of India
Branch: - Nagpur Main Branch, Nagpur
Bank Address: - Kingsway, Nagpur, 440001

Terms & Conditions

Payment Terms: - 100% in advance
Payment Mode: - Wire transfer (T/T) or PayPal, Credit/Debit Card
Packing Charges: - NIL
Delivery period: Within 7 - 10 days after placed the order
Offer validity: 15 days
Mode: Via Express Air



14.11.2022, 08:57

Почта — KOST VERONIKA [PHD0600114] — Outlook

New quotation for 25ml 100ml PPL lined hydrothermal synthesis autoclave offered by TEFIC

tefic <info@teficbio.com>

12.11.2022, Сб, 6:05

Кому: KOST VERONIKA [PHD0600114] <VERONIKA.KOST@phd.units.it>

 Вложений: 1 (683 КБ)

New quotation for 25ml 100ml PPL lined hydrothermal synthesis autoclave offered by TEFIC.pdf;

Dear Veronika,

Good day.

* 2×25ml PPL lined hydrothermal reactor price is $52 \times 2 = 104$ USD.* 2×100ml PPL lined hydrothermal reactor price is $106 \times 2 = 212$ USD.

* Extra 25ml PPL liner is 30 USD/unit.

* Extra 100ml PPL liner is 64 USD/unit.

* Gross weight is 8kg, the cheapest door to door service FEDEX express cost to Italy is 95 USD.

* Total CFR price is $104 + 212 + 30 + 64 + 95 = 505$ USD.**Please find attached our new quotation.**

Have goods in stock.

Any questions, welcome for discussion at any time.

Warmest Regards

Linda Zhu

=====

TEFIC BIOTECH CO., LIMITED

Tel: 0086-29-8905 9253

Mobile: 0086 189 0924 0926

WhatsApp: +86189 0924 0926

14.11.2022, 08:57

Почта — KOST VERONIKA [PHD0600114] — Outlook

We Chat: 189 0924 0926

Website: www.teficbiotech.com

E-mail: info@teficbio.com

=====

Kind Remind:

- 1) The payment term of TEFIC: company to company
- 2) TEFIC signature must be shown on proforma invoice

发件人: "KOST VERONIKA [PHD0600114]" <VERONIKA.KOST@phd.units.it>

发送日期: 2022-11-11 19:05:09

收件人: tefic <info@teficbio.com>

主题: RE: 25ml 100ml PPL lined hydrothermal synthesis autoclave offered by TEFIC

Dear Sir/Madam,

Good morning!

We would like to expand the order to this total:

1. 25ml PPL lined hydrothermal synthesis reactor (PPL-25 model) **2 unit**
2. 100ml PPL lined hydrothermal synthesis reactor (PPL-100 model) **2 unit**
3. 25 ml PPL liner (additional) **1 unit**
4. 100 ml PPL liner (additional) **1 unit**

May we ask for a new quotation for these items?

Thank you.

With best regards,

Veronika Kost

От: tefic <info@teficbio.com>

Отправлено: 10 ноября 2022 г. 10:37

Кому: KOST VERONIKA [PHD0600114] <VERONIKA.KOST@phd.units.it>

Тема: 25ml 100ml PPL lined hydrothermal synthesis autoclave offered by TEFIC

Dear Veronika,

Good day.

* FOB price for 25ml PPL lined hydrothermal reactor is 52 USD / whole set, include 25ml PPL liner and all stainless steel parts.

* FOB price for 100ml PPL lined hydrothermal reactor is 106 USD / whole set, include 100ml PPL liner and all stainless steel parts.

* Gross weight is 3.75kg, the cheapest door to door service FEDEX express cost to Italy is 60 USD.

* Total CFR price is $52+106+60=218$ USD.

Have goods in stock.

Any questions, welcome for discussion at any time.

Warmest Regards

Linda Zhu

=====

TEFIC BIOTECH CO., LIMITED

Tel: 0086-29-8905 9253

Mobile: 0086 189 0924 0926

WhatsApp: +86189 0924 0926

We Chat: 189 0924 0926

Website: www.teficbiotech.com

E-mail: info@teficbio.com

=====

Kind Remind:

- 1) The payment term of TEFIC: company to company
- 2) TEFIC signature must be shown on proforma invoice

发件人: "KOST VERONIKA [PHD0600114]" <VERONIKA.KOST@phd.units.it>

发送日期: 2022-11-09 22:06:40

收件人: tefic <info@teficbio.com>

主题: RE: Re:RE: 25ml 100ml PPL lined hydrothermal synthesis autoclave offered by TEFIC

Thank you for your prompt answer!

Answering your questions:

1. How many units do you need for each one?

1. 25ml PPL lined hydrothermal synthesis reactor (PPL-25 model) **1 unit**
2. 100ml PPL lined hydrothermal synthesis reactor (PPL-100 model) **1 unit**

2. Could you please provide your detailed delivery address, contact person's full name and contact phone number? We need to show these information on the new quotation, and need your delivery address to confirm the cheapest shipping cost.

The shipping address is: **Via Licio Giorgieri, 1, 34127 Trieste TS, Italy**

Looking forward to your answer.

With best regards,

Veronika

От: tefic <info@teficbio.com>

Отправлено: 9 ноября 2022 г. 14:55

Кому: KOST VERONIKA [PHD0600114] <VERONIKA.KOST@phd.units.it>

Тема: Re:RE: 25ml 100ml PPL lined hydrothermal synthesis autoclave offered by TEFIC

Dear Veronika,

Thank you so much for your kind response.

We have 2 questions should confirm with you:

1. How many units do you need for each one?

2. Could you please provide your detailed delivery address, contact person's full name and contact phone number? We need to show these information on the new quotation, and need your delivery

address to confirm the cheapest shipping cost.

Expecting your kind response as soon as possible.

Thanks.

Warmest Regards

Linda Zhu

=====

TEFIC BIOTECH CO., LIMITED

Tel: 0086-29-8905 9253

Mobile: 0086 189 0924 0926

WhatsApp: +86189 0924 0926

We Chat: 189 0924 0926

Website: www.teficbiotech.com

E-mail: info@teficbio.com

=====

Kind Remind:

- 1) The payment term of TEFIC: company to company
- 2) TEFIC signature must be shown on proforma invoice

发件人: "KOST VERONIKA [PHD0600114]" <VERONIKA.KOST@phd.units.it>

发送日期: 2022-11-09 21:46:35

收件人: tefic <info@teficbio.com>

主题: RE: 25ml 100ml PPL lined hydrothermal synthesis autoclave offered by TEFIC

Dear Sir/Madam,

Thank you for your letter.

In order to provide our office with the more detailed information, may I ask you to send the quotation including only the requested items together with the shipping cost, bank fees etc.

- 1. 25ml PPL lined hydrothermal synthesis reactor (PPL-25 model)
- 2. 100ml PPL lined hydrothermal synthesis reactor (PPL-100 model)

With best regards,

Veronika

От: tefic <info@teficbio.com>

Отправлено: 9 ноября 2022 г. 10:53

Кому: KOST VERONIKA [PHD0600114] <VERONIKA.KOST@phd.units.it>

Тема: 25ml 100ml PPL lined hydrothermal synthesis autoclave offered by TEFIC

Dear Veronika Kost,

*** FOB price for 25ml PPL lined hydrothermal reactor is 52 USD / whole set, include 25ml PPL liner and all stainless steel parts.**

*** FOB price for 100ml PPL lined hydrothermal reactor is 106 USD / whole set, include 100ml PPL liner and all stainless steel parts.**

* MOQ is 1 unit.

* Payment terms is T/T or Western union 100% in advance, also can accept online credit card payment way

* Have goods in stock, delivery date is 1~2 working days after packing

More details, please find attached our 25ml to 100ml hydrothermal synthesis autoclave official quotation.

1) TEFIC laboratory hydrothermal synthesis reactor is made of thick high quality SS304. Kettle body and kettle lid sealed with durable, reliable line sealing structure, long-term stability of the sealing effect without leakage

2) Inner uses complete PPL material, excellent acid resistance and alkali resistance.

Could you please tell us how many units do you need? so that we can confirm the cheapest door to door express cost for you as soon as possible.

Expecting your kind response.

Warmest Regards

Linda Zhu

=====

TEFIC BIOTECH CO., LIMITED

Tel: 0086-29-8905 9253

Mobile: 0086 189 0924 0926

WhatsApp: +86189 0924 0926

We Chat: 189 0924 0926

Website: www.teficbiotech.com

E-mail: info@teficbio.com

=====

Kind Remind:

- 1) The payment term of TEFIC: company to company
- 2) TEFIC signature must be shown on proforma invoice

发件人: tefic <info@teficbio.com>

发送日期: 2022-11-09 17:43:48

收件人: "KOST VERONIKA [PHD0600114]" <VERONIKA.KOST@phd.units.it>

主题: Re:Quotation_PPL lined hydrothermal synthesis autoclave

Dear Veronika Kost,

Thanks a lot for your kind inquiry on our website, this is Linda from TEFIC, professional China hydrothermal synthesis reactors manufacturer.

Our company mainly supply:

- 1) 25ml~2000ml teflon lined hydrothermal synthesis reactors, recommend operating temperature: ≤ 230 °C
- 2) 25ml~500ml PPL lined hydrothermal synthesis reactors, recommend operating temperature: ≤ 280 °C
- 3) Extra replacement teflon and PPL liners for above 2 kind of reactors

About the 25ml and 100ml PPL lined hydrothermal autoclave you inquired, will send our formal offer to you in the next E-mail.

Warmest Regards

Linda Zhu

=====

TEFIC BIOTECH CO., LIMITED

Tel: 0086-29-8905 9253

Mobile: 0086 189 0924 0926

WhatsApp: +86189 0924 0926

We Chat: 189 0924 0926

Website: www.teficbiotech.com

E-mail: info@teficbio.com

=====

Kind Remind:

- 1) The payment term of TEFIC: company to company
- 2) TEFIC signature must be shown on proforma invoice

发件人: "KOST VERONIKA [PHD0600114]" <VERONIKA.KOST@phd.units.it>

发送日期: 2022-11-08 21:22:36

收件人: "info@teficbio.com" <info@teficbio.com>

主题: Quotation_PPL lined hydrothermal synthesis autoclave

Dear Sir/Madam,

Good evening!

My name is Veronika Kost and I am a researcher at the University of Trieste (Italy). I am interested in the following items: PPL lined hydrothermal synthesis autoclave 25ml and 100ml in volume.

Could I ask you to send me the actual quotation for these items?

With best regards,

Veronika Kost



ll.4

TEFIC BIOTECH CO., LIMITED

ADD: Longfubeijun, Weiyang Palace Street, Weiyang District, Xi'an, China
Email: info@teficbio.com Att: Linda Zhu Tel: 0086-29-8905 9253 Mobile: 0086-189 0924 0926

Quotation

Quotation No:TF221110LI66

Date: November 12th, 2022

Buyer	University of Trieste	Seller	TEFIC BIOTECH CO., LIMITED
Contact	Veronika Kost	Contact	Linda Zhu
E-mail	VERONIKA.KOST@phd.units.it	Email	info@teficbio.com
Delivery Address	University of Trieste Via Licio Giorgieri, 1, 34127 Trieste TS, Italy	Tel	+86-29-8905 9253
		Address	Longfubeijun, Weiyang Palace Street, Weiyang District, Xi'an, China

DESCRIPTION OF GOODS, QUANTITY, UNIT PRICE, AMOUNT

Product name	Model	Unit Price	Qty	Total Price
25ml PPL lined hydrothermal synthesis reactor	PPL-25	USD 52	2 units	USD 104
100ml PPL lined hydrothermal synthesis reactor	PPL-100	USD 106	2 units	USD 212
Extra 25 ml PPL liner	PPL-25ml	USD 30	1 unit	USD 30
Extra 100 ml PPL liner	PPL-100ml	USD 64	1 unit	USD 64
The cheapest door to door service FEDEX express cost to Italy				USD 95
Total CFR door to door				USD 505

Please note:

1. Payment terms: T/T 100% (Buyer should bears both buyer & seller bank charges)
2. Latest delivery date: within 2 working days after receiving prepayment.
3. Warranty: one year quality warranty
4. This price is valid for one month







Linda Zhu
peter zhu 朱

TEFIC USD BANK DETAILS

Payment currency **USD**
 SWIFT/BIC **CITIKHXX** (CITIKHXXXX * If 11 characters are required)
 Account number **705049691**
 Account name **TEFIC BIOTECH CO., LIMITED**
 Bank Name: **Citibank N.A., Hong Kong Branch**
 Bank Address **3 Garden Road, Central, Hong Kong**
 Country/region **Hong Kong**
 Bank code **006**
 Branch code **391** * If paying from Hong Kong banks

Product List

Product	list	Unit price
<p>25ml PPL lined hydrothermal synthesis reactor (PPL-25 model)</p> 	<ol style="list-style-type: none"> 1. Material : Shell made of stainless steel(SS304); Liner materials is special PPL 2. Capacity: 25ml 3. Operating temperature \leq 280 °C 4. Working pressure \leq 3 Mpa (surface pressure) 5. Tempe heating and cooling speed: \leq 5°C/ min 6. Whole set reactor outer diameter: 53mm 7. Whole set reactor height: 92mm 8. PPL vessel outer diameter: 33mm 9. PPL vessel height: 61mm 10. Gross weight: 1.1kg <p>Whole set price, include stainless steel shell+ 25ml PPL liner and all accessories</p>	USD 52
<p>100ml PPL lined hydrothermal synthesis reactor (PPL-100 model)</p> 	<ol style="list-style-type: none"> 1. Material : Shell made of stainless steel(SS304); Liner materials is special PPL 2. Capacity: 100ml 3. Operating temperature: \leq 280 °C 4. Working pressure \leq 3 Mpa (surface pressure) 5. Tempe heating and cooling speed: \leq 5°C/ min 6. Whole set reactor outer diameter: 76mm 7. Whole set reactor total height: 160mm 8. PPL vessel outer diameter: 49mm 9. PPL vessel height: 91.2mm 10. Gross weight: 2.65kg <p>Whole set price, include 100ml PPL liner and all stainless steel parts.</p>	USD 106

<p>25ml PPL liner (PPL-25ml model)</p> 	<ol style="list-style-type: none"> 1. Liner materials is special PPL 2. Capacity: 25ml 3. Operating temperature: $\leq 280\text{ }^{\circ}\text{C}$ 4. Working pressure $\leq 3\text{ Mpa}$ (surface pressure) 5. Tempe heating and cooling speed: $\leq 5^{\circ}\text{C}/\text{min}$ 6. PPL vessel outer diameter: 33mm 7. PPLvessel height: 61mm 	<p>USD 30</p>
<p>100ml PPL liner (PPL-100 model)</p> 	<ol style="list-style-type: none"> 1. Material : PPL 2. Capacity: 100ml 3. Operating temperature: $\leq 280\text{ }^{\circ}\text{C}$ 4. Working pressure $\leq 3\text{ Mpa}$ (surface pressure) 5. Tempe heating and cooling speed: $\leq 5^{\circ}\text{C}/\text{min}$ 6. PPL vessel outer diameter: 49mm 7. PPL vessel height: 91.2mm 	<p>USD 64</p>

Product Structure



Hydrothermal Reactor SS304



Corrosion resistance




Customizable



Heat-Resistant

DETAILL DISPLAY



01 DETAIL DISPLAY


SS304

Acid and alkali resistance;
Corrosion resistance;
Good ductility;
Aging resistance and deformation.

02 DETAIL DISPLAY

PTFE Liner

Smooth inside and outside,
small friction coefficient for
solid materials;
Non-adhesion; Pollution-free;
Non-toxic; Easy to clean.



03 DETAIL DISPLAY

Thickened Shell

High safety
Resistance to high pressure



04 DETAIL DISPLAY

Strong sealing, Prevent leakage

Circular groove seal, safe to use.



PRODUCT DISPLAY



**TEFIC supply 10ml-2500ml teflon lined hydrothermal synthesis reactor,
25ml-500ml PPL lined hydrothermal synthesis reactor, and all specifications of
teflon liners and PPL liners**

**Thanks for your kind support all the time, warm welcome to contact us at any time,
our service at 24h /7 days.**

TEFIC official website QR code



TEFIC

TEFIC BIOTECH CO., LIMITED